

# Canon




## *Aplio i900*

Intuitive. Intelligent.  
Innovative.

CV


Cardiovascular Imaging


# Aplio i900



## Getting to the heart of the matter

Meet Aplio™ i900. Engineered to help you get the information you need to make confident decisions quickly, Aplio i900 elevates cardiovascular imaging to a new level of imaging precision, diagnostic performance and productivity.








# Crystal-clear imaging, superior versatility

From the smallest to the toughest patients, Aplio's revolutionary iBeam architecture with dramatically increased processing power\* provides unprecedented imaging clarity and definition while significantly enhancing penetration.


**iBEAM**  
ARCHITECTURE

**LIVE**  
**VOLUME**  
**IMAGING**


Aplio i900 brings best-of-breed cardiac 4D imaging to everyday clinical use. Its fully sampled matrix array transducers are especially small and lightweight for better ergonomics and patient accessibility.



Better diagnostics starts here



Aplio's ultra-wideband i-series transducers cover the same bandwidth as two conventional transducers, providing superior sensitivity and resolution for both near and far field. While helping to reduce cost, this revolutionary transducer design can provide better imaging regardless of the patient condition.



Ultra-Wideband  
Volume TEE i6SVX2



Ultra-Wideband  
Volume Sector i6SVX1



Ultra-Wideband  
Sector i6SX1




# Enjoy the perfect picture

Each of Aplio's unique imaging technologies provides you with better image quality by reducing clutter, strengthening signal and improving visualization. All functions work hand in hand with other imaging modes for greater uniformity across all applications.



- Tissue Enhancement delivers images of stunning smoothness and clarity with improved endocardial border delineation, especially in difficult patients.




- Flex-M allows you to derive anatomically correct M-mode traces from live or stored
- 2D images with the same quality as in conventional M-mode.




# Clearer vistas, faster results

Aplio's volume matrix transducer enables you to acquire high-quality 4D volumes at high framerate, while the TEE probe's small insertion head offers an improved patient experience. A host of advanced analysis tools allows for in-depth assessment of cardiac anatomy and function.



Aplio's advanced Wall Motion Tracking technology provides immediate visual and quantitative access to global and regional myocardial wall motion dynamics in 2D and 3D.



The automated MVA tool provides concise anatomic and functional assessment of the mitral valve. The function's quad display offers a clear overview of different scan planes.







ADVANCED  
WALL MOTION  
TRACKING

AUTOMATED  
MITRAL VALVE  
ASSESSMENT



At a simple touch of a button, Aplio demonstrates the mitral valve as seen by the surgeon to facilitate visual assessment of the leaflets for better surgical planning.





## Designed with our users in mind

Smaller and lighter, Aplio i900 is easy to maneuver. With over 36 cm panel height adjustment, lateral slide and a fully articulating monitor arm, Aplio i900 helps you to optimally adjust the console to virtually any scanning position.



# Aplio makes your work flow

Aplio provides a host of intelligent workflow support and automation tools, helping you to achieve rapid results with consistent high quality regardless of the patient condition.

INTELLIGENT  
PANEL

50%  
LESS HARD  
KEYS

INTERACTIVE  
TOUCH SCREEN



## Access all areas


Aplio's large, tablet-style touch screen with three interactive zones allows you to quickly browse and select the desired function, while the rest of the display remains unchanged.



# Switch to auto-pilot

Aplio's context-sensitive user interface is designed to make your imaging task simpler and quicker. While automated settings can deal with routine clinical needs, you always retain control over all imaging parameters when needed.


**INTELLIGENT  
ON-SCREEN  
NAVIGATION**




The mode-sensitive on-screen navigation for the central trackball boosts your workflow and efficiency. By visually guiding you through the exam, it allows you to adapt and operate the system within a few minutes.



- Realtime QuickScan allows you to achieve greater consistency in your exams by ensuring that superb image quality is the benchmark at all times.



- Thanks to Aplio's embedded raw data functionality you can optimize, review, analyze and report your clinical data anytime with no loss of functionality.



- A range of automated measurement and analysis tools help you increase accuracy, consistency and speed of your exams.



# Go wireless to gain better access

Aplio i-series allows you to remotely operate the system from a wireless tablet. This is especially helpful during MSK and vascular exams where it can be difficult to scan a patient and reach the panel at the same time, without losing sight on the monitor.

The wireless tablet is also ideal for scanning in mobile environments, sterile situations and for infection control to protect the system from possible contamination.




**REMOTE  
ACCESS**





Aplio's raw data capability enables you to efficiently scan and acquire images on the tablet, while you can do measurements and annotations later on the system console.



# *Aplio i900*

Intuitive. Intelligent.  
Innovative.



## ***Aplio*** *i-series*



CANON MEDICAL SYSTEMS EUROPE B.V.

<https://eu.medical.canon>

©Canon Medical Systems Corporation 2018. All rights reserved.  
Design and specifications subject to change without notice.  
Model number: TUS-AI900  
MCAUS0270ECB 2018-01 CMSE/Printed in Europe

Canon Medical Systems Corporation meets internationally recognized standards for Quality Management System ISO 9001, ISO 13485.  
Canon Medical Systems Corporation meets the Environmental Management System standard ISO 14001.

Aplio is trademark of Canon Medical Systems Corporation.

Disclaimer: Some features presented in this brochure may not be commercially available on all systems shown or may require the purchase of additional options. Please contact your local Canon Medical Systems representative for details.

***Made For life***